GENETIC ENGINEERING




Dr. Manu Smriti Singh

Department of Biotechnology


Bennett University


rDNA



ALL BEGAN WITH AN ENZYME!

- Scissors → Sticky ends
- Palindromes → Radar/Reviver





RESTRICTION ENDONUCLEASES




FIGURE 19.3

One of the first genetic engineering experiments. This diagram illustrates how Cohen and Boyer inserted an amphibian gene encoding rRNA into pSC101. The plasmid contains a single site cleaved by the restriction endonuclease *Eco*RI; it also contains *tet*^{*}, a gene which confers resistance to the antibiotic tetracycline. The rRNA-encoding gene was inserted into pSC101 by cleaving the amphibian DNA and the plasmid with *Eco*RI and allowing the complementary sequences to pair.

Cut → Restriction Endonuclease
Join → DNA Ligase

Fragments of elephant and ostrich DNA cleaved by the same endonuclease can be joined to one an-other as readily as two bacterial DNA fragments.

Genetic engineering experiments consist of four stages:

Step 1: DNA cleavage


Step 2: Production of recombinant DNA

Step 3: Cloning

Step 4: Screening

Transformation: Recombinant Plasmid+ Bacteria


- Calcium Chloride
- Electroporation




SCREENING Bacterial Genome Foreign DNA Not Inserted Bacterial Cells are resistant to tetracycline and exibit a blue color Auger contains tetracycline and Mix Plasmids Tetracycline reagent that reacts with into Resistance Gene Foreign DNA inserted at X B-Galactosidase to turn blue. **Bacterial Cells** into B-Galactosidase Coding Gene No Plasmids Present in Cell Bacterial Cells will not grow in the presence of tetracycline. B-Galactosidase Coding Gene promoter nRNA → β-galactosidase - IPTG + IPTG + X-gal + X-gal Foreign DNA Inserted Bacterial Cells are resistant to tetracycline but do not exibit the blue color.

COLONY PICKING

Screening for pUC8 recombinants




Recombinants are screened by plating onto agar containing X-gal and IPTG.





Analysis of Recombinant DNA Technology Results

Selection for Recombinant



☐ Agar dish after rec DNA Technology shows **both blue and white** (clear) colonies.


WHY both blue and white (clear) colonies?

Which of the colonies is the color we want?

If you forgot to add antibiotics in the plate, What results would you expect?




Detection Techniques





VISUALIZATION





Mem-


brane:

DNA or RNA to specific sequence (restriction fragment) DNA or RNA to specific sequence (transcript)

Use Secondary antibody to detect/amplify primary


Detect labeled probe on membrane






Sample contains specific DNA restriction fragment

Can measure fragment size and amount (single vs. repeated)



Sample contains specific RNA transcript (e.g. mRNA)

Can measure fragment size and amount (level of expression)



Sample contains specific Polypeptide

Can measure polypeptide size and amount (level of expression)



POLYMERASE CHAIN REACTION (PCR)


What is PCR?:


Use of DNA polymerase to selectively amplify a segment of DNA from a much larger sample.

Xeroxing DNA, start with one page and get many.


Examples of what PCR is used for:

- →Forensics, DNA typing from very small samples
- →Clinical diagnostics e.g. detection of HIV, detection of some microbial infections, detection of whether an individual carries a mutation predisposing them to some sort of cancer or genetic disease.
- →Research, mapping and sequencing of genomes, cloning, basic research






REAL TIME PCR (RT-PCR)



Coronavirus Testing

"Gold Standard" RT-PCR Test




>>>>


The virus is collected via a saliva sample or nose or throat swab.

Researchers break the virus, releasing its genetic material —RNA An enzyme called reverse transcriptase (RT) converts viral RNA to doublestranded DNA. Researchers use a technique called polymerase chain reaction (PCR) to create millions of copies of the converted viral DNA, making it easier to detect.


A fluorescent probe glows when it spots telltale virus DNA. This signifies a positive test.



Results














ANTIBODY DETECTION










Applications of Genetic Engineering

- Molecular biology
- Genetic disorder
- Gene therapy
- DNA fingerprinting
- Vaccines
- Pharmaceutical products

Human Insulin Production


1951: 10,000 POUNDS OF PIG PANCREASES MAKE 1 POUND OF INSULIN

TODAY:
GENETICALLY
ENGINEERED
BACTERIA PRODUCE
ANIMAL-FREE INSULIN




SUPERBUG TO CLEAR OIL SPILLS


- Super bug was developed by Anand Chakrabarty et al. in 1979.
- It is used to treat oil spills as a measure to control oil pollution.
- Petroleum products contain cycloalkenes(octane), napthenes, xylene, tolune and aromatic hydrocarbons. Since these compounds are not easily biodegradable, oil wastes become a major pollutant on the soil and water.
- Chakrabarty et al. took attempts to degrade oil wastes using micro organisms.
- They developed superbug to control oil pollution.



Anand Chakrabarty



VACCINES



FARMING




FIGURE 19.21

Transgenic rice. Developed by Swiss bioengineer Ingo Potrykus, transgenic rice offers the promise of improving the diets of people in rice-consuming developing countries, where iron and vitamin A deficiencies are a serious problem.

